

Programming an Input - Phone

Input Programming - Input 1

Input | Alarm | Phone | Email | Pre-Alarm | Advanced

Telephone Alarm: Never

Day/Night Schedule: Always Use Same List Add

Day Phone List: Phone List 1 Add

Night Phone List: Phone List 1 Add

Cancel OK

In the event of an alarm, and a voice dialer is present, there are three telephone options that can be chosen: never telephone, telephone after delay, telephone immediately.

This allows the user to set the Node so that it uses either the same or different phone lists for day and night or weekdays and weekends.

These options allow the selection of the phone list an input will call depending on the type of schedule selected.

Once a schedule has been created with its ON and OFF cycles, select normal (Day) and alternate (Night) phone lists. The default (Always Use Same List) uses only one list and does not allow the selection of a second phone list. If "Never" is selected as the telephoning option, the scheduling functions will be disabled.