

Programming an Input - Overview

The input name reflects the previously entered moniker for the given input.

Input Programming - Input 1

Input | Alarm | Phone | Email | Pre-Alarm | Advanced

Input Name: Input 1 Add Chng Del

Input Number: 1

Reads: Unused

Map Icon: Generic

Cancel OK

Add an input to the system.

Delete the selected input from the system.

Change the existing name for the given input.

In order to set up a given input you must correctly identify and select how that input is identified by the system. By clicking on Reads you can set the type of probe that is being read, the value of measurement for that probe, the scalar parameters for a linear transducer, and the ability to make an input a type of switch contact.