



MONITORING



PROTECT VALUABLE ASSETS



COMPLIANCE



NOTIFICATION



GRAPH OF EVENTS

Continuous Temperature Monitoring

Enhanced Encryption and Security Upgrade

As part of our continual goal to keep your system and software safe and reliable, Rees Scientific has integrated additional protection and security. Make your software, modules and web more secure with our latest software upgrade.



Apache

This upgrade installs the latest Apache security patch. It implements:


- Hypertext Transfer Protocol (HTTPS)
- Disables weak cipher suites
- Ensures privacy and reduces vulnerability to ransomware and other viruses



This upgrade installs and uses a new database client capable of a more secure TLS 1.2 for communications.

This upgrade:

- Reduces danger of ransomware viruses
- Protects your data privacy
- Allows you to upgrade the computer network to TLS 1.2 (to meet modern security standards)



SNMP

The software upgrade turns off the SNMP server in all of our V2 and Z3 net devices. Disabling the SNMP server in the net devices reduces vulnerability.

Upgrade at your next validation!

*Build 1021

